

Pipette Performance Data Sheet



Eppendorf Research® 3 neo Pipettes

Are you searching for a mechanical pipette that could be the best fit for you? Look no further than the new Eppendorf Research® 3 neo! Experience an efficiency boost with faster volume selection or precise one-handed tuning with the 'easy' volume adjustment setting. Improve your reproducibility with the new volume lock and temporary adjustment to increase accuracy in different workflows. Enjoy comfortable pipetting with a short pipetting button and low operating forces. And you can of course always rely on the Research 3 neo for precise and consistent pipetting, crafted with the same excellence in liquid handling that has made Eppendorf pipettes renowned worldwide.

Currently available as single-channel pipettes with nominal volumes ranging from 2.5 µL to 10 mL.

Features

- > Mechanical air-cushion pipette with 2-button operation (one operating element for volume changes / pipetting, and one for tip ejection), and customizable settings for easy, ergonomic, and reliable pipetting
- > Shorter control button with integrated volume adjustment enhances ergonomics for all and provides added comfort for small hands
- > Unique volume gear shift to either speed up volume changes by up to 40 % (fast setting) or fine-tune easily with just one hand (easy setting)
- > Volume lock to prevent unintentional volume changes
- > Pipetting performance beyond ISO 8655 requirements when used with epT.I.P.S.® pipette tips
- > Ergonomic design according to the Eppendorf PhysioCare Concept®: well-balanced, organic, lightweight, and featuring low pipetting forces
- > Spring-loaded tip cones reduce tip attachment and ejection forces, minimizing the risk of strain injuries
- > Temporary adjustment to increase accuracy for various liquid types, tip geometries, or altitude without permanently altering the factory calibration
- > Removable labeling and marking options with ColorTag pipette marking rings¹⁾
- > Full autoclavability, high UV and chemical resistance, along with easy disassembly and reassembly, facilitate versatile use and effortless decontamination and cleaning
- > Sustainability independently validated via ACT® label by NGO My Green Lab

¹⁾ sold separately

²⁾ not recommended

Applications

- > Forward pipetting
- > Reverse pipetting
- > Removal of supernatants
- > Sample mixing
- > Phase extraction
- > Filling of plates, gels, and reaction vessels

Technical Specifications

| | |
|--|---|
| Intended use | Ergonomic micropipettes for general lab use. Intended for liquid transfer within specified volumes using compatible tips. |
| Pipette type | Air-cushion system |
| Channel options | 1 |
| Volume range | 0.1 µL – 10,000 µL (11 models) |
| Volume increments | 0.002 µL – 10 µL (11 models) |
| Operation | Mechanical, 4-digit volume display |
| Volume adjustment | Unique volume gear shift with 2 modes: <ul style="list-style-type: none"> > Fast mode: Up to 40 % quicker > Easy mode: Effortless fine-tuning |
| Volume lock | Yes |
| Customization options | <ul style="list-style-type: none"> > Two volume adjustment speeds > Temp. adjustment for liquids, tips, altitude, or reverse pipetting (separate from factory adjustment) > Mark your name, lab, application, calibration date with new ColorTag marking rings¹⁾ |
| Autoclavable | Fully autoclavable (121 °C, 1 bar, 20 min) |
| UV resistance | Yes (254 nm, 30 W, 60 cm distance) |
| H ₂ O ₂ resistance | Yes |
| Consumables | <ul style="list-style-type: none"> > epT.I.P.S.® > ep Dualfilter T.I.P.S.® / ep Dualfilter T.I.P.S.® SealMax® > epT.I.P.S.® BioBased |
| Open system (compatible with 3rd party tips) | Yes ²⁾ |
| Conforms with | ISO 8655 |
| Storage accessories (sold separately) | <ul style="list-style-type: none"> > Pipette Carousel for 6 pipettes > Pipette Stand (with separate holder) > Pipette Holder for Carousel, Stand or wall mounting |
| Operating conditions | <ul style="list-style-type: none"> > Temperature: 5 °C – 40 °C > Humidity: 10 % – 95 % |
| Others | <ul style="list-style-type: none"> > Spring-loaded tip cones for easy tip attachment & ejection > High chemical resistance suitable for a wide variety of applications > Life-time tested for min. 5 years of use, incl. disassembly & reassembly, adjustment, drop tests, autoclaving > Easy maintenance, spare parts, accessories and calibration services available from Eppendorf |

Sustainability

| | |
|---------------|--|
| Certification | ACT label (pipettes), ISO 14001 (manufactory) |
| Manufactory | Assembled with 100 % renewable energy and under EU regulations in Hamburg, Germany. |
| Packaging | Fully recyclable: Made from 95% recycled cardboard, with reduced size and ink use compared to previous pipette packs. No composite materials, plastic only used for spare filters and pipette tip wrappings. |
| Maintenance | Easy disassembly and reassembly, lifetime-tested for at least 5 years of use. |
| Repairability | Reliable repair services provided by Eppendorf, and supply of spare parts for 7 years after sales stop. |

Dimensions and Storage & Transport Conditions

| | |
|--------------------------------|---|
| Length | > 225 mm (1-channel, 2.5 µL) > 235 mm (1-channel, 1,000 µL) > 210 mm (1-channel, 10 mL) |
| Weight | > 74 g (1-channel, 2.5 µL) > 80 g (1-channel, 1,000 µL) > 112 g (1-channel, 10 mL) |
| Packaging (W x H x D) | 27.5 cm x 17 cm x 9.4 cm (max., depending on model) |
| Storage & transport conditions | > Temperature: -5 °C – 45 °C > Humidity: 10 % – 95 % |
| Shipping weight | 183 g – 617 g (depending on model) |

Technical specifications are subject to change. Errors and omissions excepted.

Measurement Deviations

| Order no. | Volume range | Color code | Increment | Volume | Systematic error ¹⁾ | | Random error ¹⁾ | |
|---|--------------|---------------------------------------|-----------|----------|--------------------------------|------------|----------------------------|------------|
| Eppendorf Research® 3 neo, 1-channel, incl. epT.I.P.S.® Box 2.0 (x 96 tips) and ep Dualfilter T.I.P.S.® 'PCR clean & sterile' Rack (x 96 tips) | | | | | | | | |
| 3174 000 001 | 0.1 - 2.5 µL | ■ dark gray (for epT.I.P.S.® 10 µL) | 0.002 µL | 0.10 µL | ± 24 % | ± 0.024 µL | ± 10 % | ± 0.01 µL |
| | | | | 0.25 µL | ± 12 % | ± 0.03 µL | ± 6 % | ± 0.015 µL |
| | | | | 1.25 µL | ± 2.5 % | ± 0.031 µL | ± 1.5 % | ± 0.019 µL |
| | | | | 2.50 µL | ± 1.4 % | ± 0.035 µL | ± 0.7 % | ± 0.018 µL |
| Eppendorf Research® 3 neo, 1-channel, incl. epT.I.P.S.® Box 2.0 (x 96 tips) and ep Dualfilter T.I.P.S.® BioBased 'PCR clean & sterile' Reload (x 96 tips) | | | | | | | | |
| 3174 000 002 | 0.5 - 10 µL | ■ medium gray (for epT.I.P.S.® 20 µL) | 0.01 µL | 0.5 µL | ± 8 % | ± 0.04 µL | ± 5 % | ± 0.025 µL |
| | | | | 1 µL | ± 2.5 % | ± 0.025 µL | ± 1.8 % | ± 0.018 µL |
| | | | | 5 µL | ± 1.5 % | ± 0.075 µL | ± 0.8 % | ± 0.04 µL |
| | | | | 10 µL | ± 1 % | ± 0.1 µL | ± 0.4 % | ± 0.04 µL |
| 3174 000 003 | 1 - 20 µL | ■ light gray (for epT.I.P.S.® 20 µL) | 0.02 µL | 1 µL | ± 10 % | ± 0.1 µL | ± 3 % | ± 0.03 µL |
| | | | | 2 µL | ± 5 % | ± 0.1 µL | ± 1.5 % | ± 0.03 µL |
| | | | | 10 µL | ± 1.2 % | ± 0.12 µL | ± 0.6 % | ± 0.06 µL |
| | | | | 20 µL | ± 1 % | ± 0.2 µL | ± 0.3 % | ± 0.06 µL |
| 3174 000 004 | 1- 20 µL | | 0.02 µL | 1 µL | ± 10 % | ± 0.1 µL | ± 3 % | ± 0.03 µL |
| | | | | 2 µL | ± 5 % | ± 0.1 µL | ± 1.5 % | ± 0.03 µL |
| | | | | 10 µL | ± 1.2 % | ± 0.12 µL | ± 0.6 % | ± 0.06 µL |
| | | | | 20 µL | ± 1 % | ± 0.2 µL | ± 0.3 % | ± 0.06 µL |
| 3174 000 005 | 5 - 100 µL | ■ yellow (for epT.I.P.S.® 200 µL) | 0.1 µL | 5 µL | ± 6 % | ± 0.3 µL | ± 2 % | ± 0.1 µL |
| | | | | 10 µL | ± 3 % | ± 0.3 µL | ± 1 % | ± 0.1 µL |
| | | | | 50 µL | ± 1 % | ± 0.5 µL | ± 0.3 % | ± 0.15 µL |
| | | | | 100 µL | ± 0.8 % | ± 0.8 µL | ± 0.2 % | ± 0.2 µL |
| 3174 000 006 | 10 - 200 µL | | 0.2 µL | 10 µL | ± 5 % | ± 0.5 µL | ± 1.4 % | ± 0.14 µL |
| | | | | 20 µL | ± 2.5 % | ± 0.5 µL | ± 0.7 % | ± 0.14 µL |
| | | | | 100 µL | ± 1 % | ± 1 µL | ± 0.3 % | ± 0.3 µL |
| | | | | 200 µL | ± 0.6 % | ± 1.2 µL | ± 0.2 % | ± 0.4 µL |
| 3174 000 007 | 15 - 300 µL | ■ orange (for epT.I.P.S.® 300 µL) | 0.2 µL | 15 µL | ± 5 % | ± 0.75 µL | ± 1.4 % | ± 0.21 µL |
| | | | | 30 µL | ± 2.5 % | ± 0.75 µL | ± 0.7 % | ± 0.21 µL |
| | | | | 150 µL | ± 1 % | ± 1.5 µL | ± 0.3 % | ± 0.45 µL |
| | | | | 300 µL | ± 0.6 % | ± 1.8 µL | ± 0.2 % | ± 0.6 µL |
| 3174 000 008 | 50 - 1000 µL | ■ blue (for epT.I.P.S.® 1,000 µL) | 1 µL | 50 µL | ± 6 % | ± 3 µL | ± 1.2 % | ± 0.6 µL |
| | | | | 100 µL | ± 3 % | ± 3 µL | ± 0.6 % | ± 0.6 µL |
| | | | | 500 µL | ± 1 % | ± 5 µL | ± 0.2 % | ± 1 µL |
| | | | | 1000 µL | ± 0.6 % | ± 6 µL | ± 0.2 % | ± 2 µL |
| Eppendorf Research® 3 neo, 1-channel, incl. epT.I.P.S.® sample bag (x 10 tips) | | | | | | | | |
| 3174 000 009 | 0.1 - 2 mL | ■ red (for epT.I.P.S.® 2.5 mL) | 2 µL | 100 µL | ± 5 % | ± 5 µL | ± 1.4 % | ± 1.4 µL |
| | | | | 200 µL | ± 3 % | ± 6 µL | ± 1.2 % | ± 2.4 µL |
| | | | | 1000 µL | ± 0.8 % | ± 8 µL | ± 0.2 % | ± 2 µL |
| | | | | 2000 µL | ± 0.5 % | ± 10 µL | ± 0.2 % | ± 4 µL |
| 3174 000 010 | 0.25 - 5 mL | ■ violet (for epT.I.P.S.® 5 mL) | 5 µL | 250 µL | ± 4.8 % | ± 12 µL | ± 1.2 % | ± 3 µL |
| | | | | 500 µL | ± 2.4 % | ± 12 µL | ± 0.6 % | ± 3 µL |
| | | | | 2500 µL | ± 0.8 % | ± 20 µL | ± 0.25 % | ± 6.25 µL |
| | | | | 5000 µL | ± 0.6 % | ± 30 µL | ± 0.15 % | ± 7.5 µL |
| 3174 000 011 | 0.5 - 10 mL | ■ turquoise (for epT.I.P.S.® 10 mL) | 10 µL | 500 µL | ± 6 % | ± 30 µL | ± 1.2 % | ± 6 µL |
| | | | | 1000 µL | ± 3 % | ± 30 µL | ± 0.6 % | ± 6 µL |
| | | | | 5000 µL | ± 0.8 % | ± 40 µL | ± 0.2 % | ± 10 µL |
| | | | | 10000 µL | ± 0.6 % | ± 60 µL | ± 0.15 % | ± 15 µL |

¹⁾ The error data, according to EN ISO 8655, only apply if original Eppendorf tips are used. Technical specifications are subject to change. Errors and omissions excepted.