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Safety instructions

- Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations
- CAUTION!** For research purposes only Not suitable for therapeutic or diagnostic use in human or veterinary medicine (in accordance with the IVD Directive).
- Do not use any damaged vessel e.g. with hairline cracks or splits.
- Ensure that the cover is screwed or snapped (DIS-300-S) on tightly. **Caution!** The tube mustn’t operate without the cover.
- The tube must only be attached and removed whilst the motor is stationary.
- Ensure that the tube is firmly attached to the bayonet lock connector on the drive unit prior to operating the device.
- Ensure that the tube is able to withstand the chemical properties of the selected medium by filling the tube and storing for a suitable period time.
- Always open the tube carefully after use as the media in the tube may heat up due to transfer of energy during operation, leading to pressurization of the container: risk of material spraying, protective equipment must be worn.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.
- The temperature of the material must not exceed 40 °C.
- Note the operating instructions of the disperser unit.

Product information

All tubes consist of a mixing chamber and agitator, plus depending on the variant an additional pulverising tool.

Nomenclature	ST	BMT	DT	DIS
Tool	Stirring unit	Stirring unit, addition-ally with balls	dispersing unit	Dissolving unit
Volume: 20/50/300 ml <p>Materials: S: stainless steel / G: glass / T: titanium</p>	ST-20 <p>ST-50</p>	BMT-20-S <p>BMT-20-G</p> <p>BMT-50-S</p> <p>BMT-50-G</p>	DT-20 <p>DT-20 eco</p> <p>DT-50</p> <p>DT-50 eco</p>	DIS-50-S-M <p>DIS-50-T-M</p> <p>DIS-300-S-M</p> <p>DIS-300-S</p>

The tubes are ideal for applications where no cross-contamination is permissible.

The tubes are intended exclusively for disposable use. Re-use can seriously damage both the tube and the drive.

The tubes are not free of traces of DNA or RNA.

Application instruction

Store the tubes in a dry place, in their packaging and at room temperature.

Note the use-by date of the tubes.

Protect the tubes against ultraviolet light.

Don’t turn the tube after piercing the membrane in the screw cap, then the leak-tightness is ensured no longer.

Always position the sample in the centre of the dispersing unit.

The solid bodies in samples must not be frozen or too hard.

Dispersing causes the medium to heat up.

Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute.

Longer application times bring no improvements, increase however the sample temperature substantially. In addition, excessive operating time can lead to material fatigue, which can cause the tube to leak.

Dispose of the tubes after use according to the regulations. Observe the laboratory regulations.

	ST	BMT	DT	DIS
Mixing chamber and cover	Polypropylene (PP), DIS-300-S cover: Polyethylen (PE)			
Membrane	Polypropylene (PP)			
Tumbling bar	Polyetheretherketone (PEEK) (DIS 300: PEEK+PTFE+graphite +carbon fibre)			
Balls	-	S: Stainless steel AISI304 <p>G: Soda lime glass</p>	-	-
Dispersing unit	-	-	Polyetheretherketone with polytetrafluoroethylene (PEEK+PTFE) <p>eco: Polyethersulfone with polytetrafluoroethylene (PES+PTFE)</p>	-
Dissolver	-	-		S: Stainless steel AISI304 <p>T: Titanium</p>
FDA conformity	Yes			
ICH Q3D / USP chapter 232 / 233 conformity	No	No	No	DIS-50-T-M
Dry operation	Yes	Yes	No	No
Sterile variants	ST-20-M-gamma <p>ST-50-M-gamma</p>	BMT-20-S-M-gamma <p>BMT-50-S-M-gamma</p>	DT-20-M-gamma <p>DT-50-M-gamma</p>	-
Chemical resistance	As for the listed materials			

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