

# Combination of epT.I.P.S.<sup>®</sup> Pipette Tips with Eppendorf Research<sup>®</sup> plus Pipettes

Eppendorf Research <sup>®</sup> plus/ Eppendorf Research <sup>®</sup> plus Move It <sup>®</sup>	0.1 - 10 µL, dark gray	0.1 - 10 µL M, medium gray	0.1 - 20 µL, medium gray	0.5 - 20 µL L, light gray	0.1 - 20 µL, light pink	2 - 20 µL, yellow	2 - 100 µL, yellow	5 - 100 µL, light yellow	2 - 200 µL, yellow	20 - 300 µL, orange	50 - 1,000 µL, blue	50 - 1,250 µL, green	50 - 1,250 µL L, dark green	0.1 - 2.5 mL, red	0.1 - 5 mL, violet	0.2 - 5 mL L, violet	0.5 - 10 mL, turquoise	0.5 - 10 mL L, turquoise
epT.I.P.S. <sup>®</sup> BioBased	■		■	■					■	■	■	■	■	■	■	■		
epT.I.P.S. <sup>®</sup> Standard	■		■	■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. <sup>®</sup> BioBased Reloads			■ ■			■ ■	■		■	■	■	■	■					
ep Dualfilter T.I.P.S. <sup>®</sup>	■	■		■	■	■	■	■	■	■	■	■	■	■	■	■		■
ep Dualfilter T.I.P.S. <sup>®</sup> BioBased Reloads	■	■		■		■	■		■	■	■	■	■					
ep Dualfilter T.I.P.S. <sup>®</sup> SealMax <sup>®</sup>	■			■			■		■	■	■							
ep Dualfilter T.I.P.S. <sup>®</sup> SealMax <sup>®</sup> BioBased Reloads				■			■		■	■	■							
<b>Fixed volume, 1-channel pipettes</b>																		
10 µL	●	■	■	●														
20 µL	■		●	■														
10 µL	■					■	■		■	●								
20 µL	■					■	■		■	●								
25 µL/50 µL/100 µL	■						■		■	●								
200 µL	■								■	●								
200 µL/250 µL/500 µL/1,000 µL	■										■	●	●					
<b>Variable volume, 1-, 8- and 12-channel pipettes</b>																		
0.1 - 2.5 µL	■	■	●	●														
0.5 - 10 µL	■	●	■	■	●													
2 - 20 µL	■	●	●	●	■													
2 - 20 µL						■	■		■	●								
10 - 100 µL	■					●	■		■	●								
20 - 200 µL	■					●	●		■	●								
30 - 300 µL	■					●	●		●	■								
100 - 1,000 µL	■										■	●	●					
120 - 1,200 µL	■											●	■					
0.25 - 2.5 mL	■													■				
0.5 - 5 mL	■														■	●		
1 - 10 mL	■																■	●
<b>Variable volume, 16- and 24-channel pipettes</b>																		
1 - 20 µL	■				■													
5 - 100 µL	■							■										
<b>Move It<sup>®</sup>, variable volume, 4-, 6-, 8- and 12-channel pipettes</b>																		
1 - 20 µL	■				■													
5 - 100 µL	■							■										
30 - 300 µL	■					●	●		■									
120 - 1,200 µL	■											●	■					