



**DURAN
WHEATON
KIMBLE**

Excellence in your hands

DURAN® PURE GLS 80®

DEVELOPED FOR PHARMA


DURAN® PURE is a range of high-quality bottles and closures developed for the needs of the pharmaceutical industry. **DURAN® PURE** meets the requirements of the international pharmacopeias and is available with the certificates required by the regulatory authorities.

FEATURES

- **EP/USP Type I borosilicate glass**
- **Size range 500 ml - 20 l**
- **GLS 80®** thread for fast closing & opening
- **Protective cover** prevents contamination
- **Single batch delivery** CoA available
- **100% traceability**
- **100% visual inspection**
- **Change control**
- **Drug Master File**

APPLICATION AREAS

- Manufacturing of pharmaceutical products
- Storage of bulk drug intermediates
- Transportation of APIs



DURAN® PURE GLS 80® Bottles

Article number clear glass	Article number amber coloured	Volume [ml]	Bottles per box [pcs.]	Bottles per pallet [pcs.]
81 860 44 03	81 866 44 09	500	10	280
81 860 54 08	81 866 54 05	1 000	10	210
81 860 63 01	81 866 63 07	2 000	10	120
81 860 73 06	81 866 73 03	5 000	1	48
81 860 86 02	81 866 86 08	10 000	1	30
81 860 91 01	81 866 91 07	20 000	1	12

ALSO AVAILABLE: DURAN® PURE bottles and screw caps with GL 45 thread

DURAN® PURE GLS 80® Screw Cap

The DURAN® PURE GLS 80® Screw Cap made from PSU is chemically and thermally very resistant and compatible with autoclaving and other common sterilization methods. The product includes a silicone seal which is coated with highly inert PTFE.

FEATURES

- Cap made from PSU with PTFE coated silicone seal
- GLS 80® thread for fast closing & opening
- Single batch delivery CoA available
- 100% traceability
- Change control



Article number	Screw Caps per PE bag [pcs.]	Screw Caps per pallet [pcs.]
89 116 91 09	1	1435



DWK Life Sciences GmbH
Hattenbergstraße 10
55122 Mainz, Germany
www.DWK.com

