

Bottle Compatibility Chart

Bottles					
Bottle Manufacturer	Volume (mL)	Order No.	Size (mm) Diameter x Height w. cap	Bottom shape	Material
BD®/Corning®	175	352076 ¹⁾	60.7 x 120.4	conical	PP
BD®/Corning®	225	352075 ¹⁾	61 x 139,4	conical	PP
Nalgene®	175	3143-0175 ¹⁾	61,8 x 139,4	conical	PPCO
Nalgene®	250	3144-0175 ¹⁾	61,8 x 140,5	conical	PC
Nalgene®	250	3145-0175 ¹⁾	61.5 x 135.6	conical	PS
Nalgene®	250	3120-0250	61.8 x 125.5	flat	PPCO
Nalgene®	250	3121-0250	61,2 x 131,6	flat	HDPE
Nalgene®	250	3122-0250	61,8 x 127,6	flat	PC
Nalgene®	250	3127-0250	60 x 126,5	flat	FEP
Nalgene®	250	3140-0250	61,8 x 134,4	flat	PC
Nalgene®	500	3141-0500	69,5 x 160	flat	PP
Nalgene®	1000	3120-1010	98 x 168.4	flat	PPCO + sealing cap
Nalgene®	1000	3122-1010	98.1 x 169.7	flat	PPCO
Nunc®	200	376813	59,9 x 137,2	conical	PP
Nunc®	250	376814	60.2 x 145	conical	PP
Corning®	250	430776	60 x 161,3	conical	PP
Corning®	250	431840	60 x 130	flat	PP + sealing cap
Corning®	250	431841	60 x 128	flat	PP w/o sealing
Corning®	250	431842	62 x 131	flat	PC + sealing cap
Corning®	250	431843	62 x 129	flat	PO w/o sealing
Eppendorf®	250	5721411013	62 x 124	flat	PC
Eppendorf®	250	5721411012	62 x 124	flat	PP
Eppendorf®	400	5810 729 009	80 x 133	flat	PP
Eppendorf®	400	5910 707 002	81 x 117	flat	PP
Eppendorf®	500	5820 707 003	83 x 134	flat	PP
Eppendorf®	750	5820 708 000	102 x 133	flat	PP
Eppendorf®	1000	5910 770 006	118 x 133.1	flat	PP
Corning®	500	431123	95 x 147	conical	PP
Corning®	500	431844	68 x 168	flat	PP
Corning®	500	431847	69 x 169	flat	PC
Beckman Coulter®	500	355605	69 x 160	flat	PC
Beckman Coulter®	500	355607	69 x 160	flat	PP
TPP®	600	87600	100 x 183	conical	PP

Eppendorf Rotors for Centrifuge SpinPro® 6 R		
S-2-x Universal ID	with AT Caps*	FA-6x250 ID
■	■	■
■	■	–
■	–	–
■	–	–
■	■	–
■	■	■
■	■	■
■	■	■
■	■	■
■	■	■
–	–	–
–	–	–
■	■	–
–	–	–
–	–	–
■	■	■
■	■	■
■	■	■
■	■	■
–	–	–
–	–	–
–	–	–
■	■	–
■	■	–
■	–	–
–	–	–
–	–	–
–	–	–

Technical specifications subject to change. Errors and omissions excepted.

■ Yes * Aerosol tight centrifugation with cap possible

1) Please note that these bottles, when centrifuged in a flat bottom shape adapter, need to be used together with a conical form insert provided by the manufacturer.

www.eppendorf.com/centrifugation

Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · Barkhausenweg 1 · 22339 Hamburg · Germany
 E-mail: eppendorf@eppendorf.com · www.eppendorf.com

Beckman Coulter® is a registered trademark of Beckman Coulter, Inc., USA. Corning® is a registered trademark of Corning, Inc, USA. Nalgene® is a registered trademark of Nalge Nunc International Corporation.
 Nunc® is a registered trademark of Nunc A/S Corporation, Denmark. TPP® is a registered trademark of TPP Techno Plastic Products AG, Switzerland.
 Eppendorf®, SpinPro® and the Eppendorf Brand Design are registered trademarks of Eppendorf SE, Germany · All rights reserved, including graphics and images · EN5/1125/HQ · Copyright © 2025 by Eppendorf SE.