

HTL™ Discovery™ Comfort single-channel pipettor shafts come in red, yellow, green, blue, or white depending on the volume range of the pipettor and the corresponding tip size.

HTL Discovery Comfort Single-channel pipettors

Model	Cat. No.	Color Code	Volume (µL)	A (%)	P (%)	Non-filter Tips
D2	4051-D	●	Min. 0.1 1.0 Max. 2.0	±40.0 ± 2.7 ± 1.5	≤12.0 ≤ 1.3 ≤ 0.7	10 µL
D10	4052-D	●	Min. 0.5 5.0 Max. 10.0	± 4.0 ± 1.5 ± 1.0	≤ 2.8 ≤ 0.7 ≤ 0.4	
D20	4053-D	●	Min. 2 10 Max. 20	± 3.5 ± 1.0 ± 0.8	≤ 1.5 ≤ 0.5 ≤ 0.3	
D50	4057-D	●	Min. 5 25 Max. 50	± 3.5 ± 1.0 ± 0.8	≤ 1.3 ≤ 0.5 ≤ 0.3	200 µL
D100	4054-D	●	Min. 10 50 Max. 100	± 3.0 ± 1.0 ± 0.8	≤ 1.0 ≤ 0.3 ≤ 0.2	
D200	4055-D	●	Min. 20 100 Max. 200	± 2.0 ± 1.0 ± 0.6	≤ 0.7 ≤ 0.3 ≤ 0.2	
D250	4050-D	●	Min. 50 125 Max. 250	± 1.0 ± 0.8 ± 0.6	≤ 0.4 ≤ 0.3 ≤ 0.3	300 µL
D1000	4056-D	●	Min. 100 500 Max. 1,000	± 2.5 ± 0.8 ± 0.6	≤ 0.6 ≤ 0.3 ≤ 0.2	1,000 µL
D5000	4058-D	○	Min. 500 2,500 Max. 5,000	± 1.2 ± 0.6 ± 0.5	≤ 0.5 ≤ 0.2 ≤0.15	5,000 µL
D10000	4059-D	○	Min. 1,000 5,000 Max. 10,000	± 2.5 ± 0.8 ± 0.5	≤ 0.6 ≤ 0.3 ≤ 0.2	10,000 µL

HTL Discovery Comfort Multi-channel pipettors

Model	Cat. No.	Color Code	Volume (µL)	A (%)	P (%)	Non-filter Tips
DV8-10 DV12-10	5121 5125	●	Min. 0.5 5 Max. 10	±10.0 ± 4.0 ± 2.0	≤ 8.0 ≤ 2.0 ≤ 1.2	10 µL
DV8-50 DV12-50	5122 5126	●	Min. 5 25 Max. 50	± 4.0 ± 3.0 ± 1.6	≤ 2.5 ≤ 1.2 ≤ 0.6	200 µL
DV8-200 DV12-200	5123 5127	●	Min. 20 100 Max. 200	± 3.0 ± 1.5 ± 1.0	≤ 1.5 ≤ 0.8 ≤ 0.6	
DV8-300 DV12-300	5124 5128	●	Min. 50 150 Max. 300	± 1.6 ± 1.2 ± 1.0	≤ 1.5 ≤ 1.0 ≤ 0.6	300 µL

NOTE: Actual products may vary from the images shown.

